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Tel. +34 854 53 63 74

satcliente@lrdiagnostico.com

www.lrdiagnostico.com

VACUTEST[®]
kima

142600

KIMASED RACK FOR 10 TUBES READING

The KIMASED RACK features a graduated scale designed to obtain results in 30'.
 The result is expressed in mm/hour as in the traditional Westergren system.

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142625

KIMASED AUTO 20

ESR autoanalyzer 20 tubes positions.
 Testing time 24 minutes - 50 samples/hour.
 With built-in printer.
 The result is expressed in mm/hour as in the traditional Westergren system.

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142645

KIMASED AUTO 60

ESR auto-analyzer 60 tubes positions.
 Testing time 24 minutes - 150 samples/hour.
 With built-in printer.
 The result is expressed in mm/hour as in the traditional Westergren system.

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* KIMASED is supplied in aluminium bags which preserve the shelf life for **18 months**. Once the tubes are taken from the bag must be used in 4 months (PAO). If tubes are kept in the bags **closed again** and **sealed** the shelf life is guaranteed for 18 months.

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CE **STERILE R**
 interior


PAO
 4M

Medical diagnostic in vitro device, Directive 98/79/CE

REF	14250	Vacuum	1.6 ml
LOT		Q.TY	
OP		DAY	

Recommended use
Mixing indications:

- 1 - immediately after blood collection gently invert the sample 8-10 times.
- 2 - before testing mix again manually or mechanically for 3-5 minutes.


Sample preservation: at room temperature (maximum 4 hours).

Storage temperatures	Up to 24 °C	at 2-4 °C
Maximum preservation time	4 hours	≤ 24 hours

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For further information related to ESR test see CLSI document (ex NCCLS) H02-A5.

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
Recommended use

Mixing indications:

- 1 - immediately after blood collection gently invert the sample 6-8 times.
- 2 - before testing mix again manually or mechanically for 2 minutes.

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