



HbA1c results in minutes

Informed conversations at the point of care

The A1CNow[®] system from PTS Diagnostics provides healthcare professionals with a fast and easy way of getting HbA1c results—allowing informed conversations about how patients are managing their diabetes in minutes, not days.

For more information, visit www.a1cnow.com.

Key Features

- National Glycohemoglobin Standardization Program certified
- CLIA-waived
- Tests may be reimbursable



Distribuido por:




Tel. +34 854 53 63 74
satcliente@lrdiagnostico.com
www.lrdiagnostico.com

HbA1c results in minutes

Informed conversations at the point of care



A1CNow[®]+ system
(10-count test kit)
3024



A1CNow[®]+ system
(20-count test kit)
3021

A1CNow[®]+

More efficient than the lab

- Results in 5 minutes
- Small (5µL) blood sample from fingerstick
- Portable – use in multiple exam rooms
- Easy to use – minimal training required
- Tests may be reimbursable



Website: ptsdiagnostics.com
Office: +1 (317) 870-5610 FAX: +1 (317) 870-5608
Toll-Free in the US: 1 (877) 870-5610