

Atellica VTLi Patient-side Immunoassay Analyzer

A vital leap forward in cardiac testing

Speed meets accuracy where it matters most—
the wait is over for hs-cTnI at the point of care.

siemens-healthineers.com/vtli



Atellica® VTLi Analyzer Specifications

System description		Patient-side immunoassay analyzer
System menu	Parameter	High-sensitivity cardiac troponin I (hs-cTnI) BNP*
	Measurement range	LoD—1250 ng/L [†] LoD plasma 1.2 ng/L [†] LoD whole blood 1.6 ng/L [†]
Sample types	Whole blood	Venous (lithium heparin [†]), capillary
	Plasma	Venous (lithium heparin [†])
Sample size		30–100 µL
Time to result		~8 minutes
Calibration		Reagent cartridges contain calibration parameters valid on RFID tag until expiry. Optical system calibration is automatically performed during every test
Quality control		Sero PATHONORM™ Cardiac Acute Liq
Integrated barcode scanner		For patient ID, operator ID, and bottled QC; 1D (linear) and 2D bar-code symbologies
External interfaces		POCT1A
Docking station power requirements	Input	100–240 VAC, 50–60 Hz, 660 mA
	Output	19 VDC, 1570 mA
	Classification	Class II

Distribuido por:



Tel. +34 854 53 63 74
satcliente@lrdiagnostico.com
www.lrdiagnostico.com

SIEMENS
Healthineers

Atellica VTLi Analyzer Specifications (continued)

Operation on battery		~60 tests on a fully charged battery	
Safety		IEC 61010-1, IEC 61010-2-101	
EMC		IEC 61326-1, IEC 61326-2-6	
Data security		Patient data anonymization TLS 1.2-encrypted end-to-end communication using client and server certificates Configurable level of user authentication No storage of passwords, no hard-coded passwords No user access to operating system	
Hardware	Operating system	Microsoft WEC2013	
	Processor	NXP i.MX 6SoloX	
	Memory	512 MB RAM	
WI-FI		802.11 b/g/n/ac 2.4 GHz	
WI-FI security		WPA/WPA2 Personal (pre-shared key)	
Data capacity	Patient/QA tests	100 tests	
	Operators	Not limited	
Reagent cartridge	Storage	2–8°C	
	Shelf life	7 months	
	Size of box (W x D x H)	152 x 95 x 167 mm	
	Cartridges per box	24	
System dimensions	Analyzer	Length	25 cm
		Height	5.2 cm
		Width	8.5 cm
		Weight (incl. battery)	780 g
		Weight	580 g
	Display	4.3-inch color LCD touchscreen with backlight	
	Docking station	Length	29 cm
		Height	6 cm
		Width	10 cm
		Weight	460 g
Environmental requirements	Operating temperature	5–27°C	
	Humidity	20–80% (noncondensing)	
	Air pressure	70–110 kPa	
	Storage temperature	–20 to 55°C	
	Storage humidity	5–85% (noncondensing)	

Atellica and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

Not available for sale in the U.S. Product availability will vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Legal Manufacturer

Siemens Healthcare Diagnostics Inc.
Laboratory Diagnostics
511 Benedict Avenue
Tarrytown, NY 10591-5005
USA
Phone: +1 914-631-8000

Local Contact Information

Siemens Healthcare Diagnostics Inc.
Point of Care Diagnostics
2 Edgewater Drive
Norwood, MA 02062-4637
USA
Phone: +1 781-269-3000