

# BS-2000 Modular System

Mindray clinical chemistry solution for high volume laboratories

## Enhanced Processing Capability

2000 tests/hour and 2200 tests/hour with ISE for each module

140 positions in sample carousel and 300 positions in SDM

Supports sample continuous loading and offloading

Reagent volume: 10~200µl

Minimum reaction volume: 80µl

Absorbance range: 0~3.5 A (10mm conversion)

System auto-start function



## Advanced Software

Supports touch screen interface

Multi-functional and easy to use

Auto re-run function

Auto serum index

Real-time QC status monitoring

Test summary function

Reflex function



## Dedicated System Pack Reagents

48 parameters of dedicated reagents available

More than 17 others are coming

Original calibrators with metrological traceability for all parameters

Original controls for all parameters

## Flexible Scalability

Offers scalability upgrade for laboratories as test volume increases

Flexible configurations with single and dual modules




Connectivity with Mindray's immunoassay platform



# BS-2000 Modular System

- Discrete, random access, fully automated
- Constant 2000 tests per hour, up to 2200 tests per hour with ISE
- Up to 67 on-board chemistries & 3 ions
- Flexible sample carousel and sample delivery module (SDM) configuration
- Minimum reaction volume 80 µl
- Grating optical system with 13 wavelengths
- Intelligent auto-recovery from probe collision
- Dedicated system pack reagents
- Flexible scalability within different modules



<p>Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 E-mail: intl-market@mindray.com www.mindray.com</p> <p><b>Mindray is listed on the NYSE under the symbol "MR"</b></p>	<p><b>mindray</b> is a trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. ©2013 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved. Specifications subject to changes without prior notice. P/N: ENG-B52000M-210285x2-20130303</p>	
	 	

Distribuido por:



Tel. +34 854 53 63 74

satcliente@lrdiagnostico.com

www.lrdiagnostico.com

**mindray**  
healthcare within reach