

Mechanical detection method, insensitive to abnormal samples that are hemolysis, lipemic, icteric etc.



C2000-4

C2000-2



Technical specifications

Testing method	Coagulation method (mechanical)	Accuracy	±3%
Testing items	PT, APTT, TT, FIB etc.	Testing Channel	2/4
Sample volume	≤200μL	Incubation Channel	6/16
Reagent volume	50μL-200μL	Reagent Position	3/4
Timing start	Electronic Start-Pipette	Printing	Built-in Printer for C2000-2/4
Channel consistency	4%		

Mindray Building, Keji 12th Road South,
High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R. China
Tel: +86 755 8188 8998 Fax: +86 755 26582680
E-mail: intl-market@mindray.com www.mindray.com

mindray (— — —) are registered trademarks or trademarks owned by Shenzhen Mindray Bio-medical Electronics Co., Ltd.
©2015 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved. Specifications subject to changes without prior notice.
P/N:ENG-C2000-2-4-210285X2P-20170302

mindray

Distribuido por:



Tel. +34 854 53 63 74
satcliente@lrdiagnostico.com
www.lrdiagnostico.com

C2000-2/4

Semi-Auto Coagulation Analyzer

Dual magnetic circuit beads method
—A promise of reliable result

mindray
healthcare within reach