

MW-12A

Microplate Washer

- Compatible with flat, U- or V-bottom plates
- 2 wash modes for 8- and 12-well strips
- Automatic bottom detection and collision protection
- Variable Shaking and soaking time
- Up to 100 user-defined protocols
- Sweep function with aspiration
- On-board automatic self-diagnostic program
- Automatic liquid monitoring system



Distribuido por:



Tel. +34 854 53 63 74
 satcliente@lrdiagnostico.com
 www.lrdiagnostico.com

mindray
 healthcare within reach

MW-12A

Microplate Washer

Technical Specifications

Wash programs	up to 100
Wash cycles	0 to 9
Dispensing volume	50 up to 3000 μ L
Bottom collision protection	automatic
Residual volume per well	< 2 μ L
Soaking time	0 to 300 seconds
Shaking time	0 to 300 seconds
Multiple options	1 x 8, 1x 12
Standard bottle set	2 x 2L (detergent), 1 x 2L (Waste) 1 x 2L (distilled water)
Power supply	AC 100-240V, 50/60Hz
Temperature	5°C~35°C
Humidity	15~85%
Dimension	400mm (L) x 330mm (W) x 145mm (H)
Weight	7.2 Kg



DISTRIBUTOR:



MINDRAY is a trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Specifications subject to changes without prior notice.
 © 2010 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.
 P/N: ENG-MW12A-210285x2-20100203

mindray

Mindray is listed on the NYSE under the symbol "MR"
 Mindray Building, Keji 12th Road South, High-tech Industrial Park,
 Nanshan, Shenzhen 518057, P.R. China
 Tel: +86 755 26582888 Fax: +86 755 26582680
 E-mail: intl-market@mindray.com Website: www.mindray.com